

## **Recommendations**

1. Haemoglobin sulphate (SLS) method provides a reliable routine method with the advantage of non-hazardous reagent.
2. We recommend another study using sodium lauryl sulfate in automated analysers.

## **References**

1. ICSH: Recommendations for reference method for haemoglobinometry in human blood (ICSH standard Ep 6/2: 1977) and specification for international haemoglobinocyanide reference preparation (ICSH standard Ep 6/3: 1977). J Clin Pathol, 1978; 13: 139-43.
2. Recommendation for reference method for haemoglobinometry in human blood (ICSH standard 1986) and specifications for international haemoglobincyanide reference preparation (3<sup>rd</sup> edition). Clin Lab Haematol, 1987; 9: 73-9.
3. Hoffbrand AV, Pettit JE, and Moss PAH: Essential Haematology, 4<sup>th</sup> edition, Blackwell Science, London, 2001; p.p. 14-7.
4. Sunderman FW, Fate RPM, Fadyen DAM, Steven GF, and Copeland BE: Symposium on Clinical Haemoglobinometry. Am J Clinical Pathol, 1953; 23: 519-98.
5. Sunderman FW: Status of clinical haemoglobinometry in the United States. Am J Clin Pathol, 1965; 93: 9-15.
6. Wolf HU, Lang W, and Zarder R: Alkaline haematin D-575, a new tool for the determination of haemoglobin as alternative to the cyanhaemoglobin method. II. Standardization of the method using pure chlorohaemin. Clin Chim Acta, 1984; 136: 95-104.
7. Severini G: The extinction coefficient of human hemoglobincyanide as determined by nitrogen analysis. Clin Chem Acta, 1966; 14: 376-7.
8. Hamaguchi Y: Overview of the principles of Sysmex's Haemoglobinometry. Sysmex J Internat, 1999; 7(1): 45-51.
9. Matsubara T: Evaluation of the internationally standardized method for haemoglobinometry. Clin Chem Acta, 1969; 23: 427-30.
10. Oshiro I, et al.: New manual and automatic method of haemoglobin determination by using sodium lauryl sulfate (SLS). Clin Pathol, 1981; 29: Abstract.
11. Oshiro I, Tuke Naka T, and Maeda J: New method for haemoglobin determination by using sodium lauryl sulfate (SLS). Clin Biochem, 1982; 15(2): Abstract.
12. Dacie JV, and Lewis SM: Practical Haematology, 8<sup>th</sup> edition, Churchill Living Stone, Edinburgh, 1995; p.p. 49-54.
13. Lewis SM, Garvey B, Manning R, Sharp SA, and Wardle J: Lauryl sulfate haemoglobin: a non hazardous substitute for HiCN in haemoglobinometry. Clin Lab Haemat, 1991; 13: 279-90.
14. ICSH: Recommendation for reference method for haemoglobinometry in human blood. Clin Lab Haematol, 1986; 9: 73-9.
15. Dacie JV, and Lewis SM: Practical haematology, 7<sup>th</sup> edition. Churchill Living hstane, Edinburgh, 1991; p.p. 37-41.
16. Karsan A, Maclaren I, Conn D, and Wadsworth L: An evaluation of haemoglobin determination using sodium lauryl sulfate. Am J Clin Pathol, 1993; 100(2): 123-6.
17. Fujiware C, et al.: The reagent sulfolysers for Hb measurement by haematology analyzers Sysmex J Internat, 1990; 13: Abstract.
18. Harmening D: Modern Blood banking and transfusion practice, 2<sup>nd</sup> edition, F. A. Davis Company Philadelphia, 1989; p.p. 65.