

Table 2: Parameters and Results

Wavelength(nm)	Parameters	Color of tattoo	Type of tattoo	Results
1064	Pulse width 10ns Spot size (3mm) R.R 5Hz Fluence 6.3-10.1j/cm ²	Blue	Professional	(20) Persons required (5-6) treatment sessions but all patients still have shadow.
			amateur	(38) Persons required (3-4) treatment sessions but all patients still have slight shadow.
1064	Pulse width 10ns Spot size (3mm) R.R 5Hz Fluence 6.3-10.1j/cm ²	Black	Professional	(10) persons required (3-5) treatment sessions for removing completely except (6) Patients still have cosmetically acceptable slight shadow.
			amateur	(6) Persons required (2-4) treatment sessions for removing completely tattoos except (2) patients still have cosmetically acceptable slight shadow.
1064	Pulse width 10ns Spot size (3mm) R.R 5Hz Fluence 6.3-10.1j/cm ²	Green	Professional	(8) Persons not response to treatment after (9-10).secessions.
			amateur	(3) Persons not response to treatment after (9-10).secessions except (1) patient slight faint in color after (9) sessions.
532	Pulse width 7ns Spot size (2mm) R.R 10Hz Fluence 7.3-10.3j/cm ²	Red	Professional And amateur	Always the red color tattoo in companion with other color (Blue, black, green). (14) Red color tattoos all of them professional required (2-3) treatment sessions to be removed completely.